Disposable Gloves

Model: LD01/LD09T/LDT09/LD09E/LD49

Konformitätserklärung und vollständige Dokumentation:
Declaration of conformity and complete documentation:
Déclaration de conformité et documentation complète:
Dichiarazione di conformità e documentazione completa:
Declaración de conformidad y documentación completa:
Verklaring van overeenstemming en volledige documentatie:
Deklaracja zgodności i pełna dokumentacja:
<https://www.lanonpro.com/pages/declarations-of-conformity>

PSA-Konformitätsbewertungsstelle:
PPE Conformity Assessment Body:
Organisme évaluation conformité EPI:
Organismo valutazione conformità DPI:
Organismo Evaluación Conformidad EPI:
PBM-conformiteitsbeoordelingsinstantie:
Jednostka oceniająca zgodność i środków ochrony
indywidualnej: **Notified body:**
SATRA Technology Europe Ltd (2777)
Bracetown Business Park, Clonee
Dublin 15015 YN2P Ireland
Certificate: 2777/12152-03/E10-01

DE

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EN

These are non-sterile protective gloves, intended for use when there is a need for protection against dangerous substances and infective agents. **These products are classed as Category III Personal Protective Equipment (PPE) by the PPE Regulation (EU) 2016/425 and PPE Regulation (EU) 2016/425 on personal protective equipment, as amended to apply in GB.**

They have been shown to comply with this regulation through the Harmonised / Designated Standard(s): EN ISO 21420:2020, EN ISO 374-1:2016+A1:2018, EN ISO 374-5:2016 and EN ISO 374-4:2019. Degradation results indicate the change in puncture resistance of the gloves after exposure to the challenge chemicals EN ISO 374-5:2016. The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen.

Use and Warning: This information does not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals. The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only and relates only to the chemical tested. It can be different if the chemical is used in a mixture. It is recommended to check that the gloves are suitable for the intended use because the conditions at the work place may differ from the type test depending on temperature, abrasion and degradation. When used, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves. Before usage, inspect the gloves for any defect of imperfections.

Donning: Select the appropriate size glove & don't use damaged gloves or imperfections. Remove one glove from the box by the cuff and place on the opposite hand. Remove a further glove from the box with the hand wearing a glove and put onto the other hand.

Doffing: Slide one glove off down to mid-palm from the cuff, repeat with the second glove turning the glove inside out. Hold both gloves without touching the contaminated wall & dispose them according to the environment requirement. The gloves have no mechanical protection. The product complies with FDA 21C.F.R Safe for food use.

The EU Declaration of Conformity could be found on website:
<https://www.lanonpro.com/pages/declarations-of-conformity>

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LANON is a registered trademark of LANON Protection LLC.
LANON Protection Technology Co., Ltd
LANON Protection LLC